

- Tips and Tricks -

1.

Application without pinching

To avoid pinching, place your index finger between the strap / buckle and the skin when tightening.



2.

Loosen first, then slide open

To relax the pressure and loosen the stowage, simply tilt the "lid" of the daisygrip towards the band.

To open:
Slide buckle to the side with your thumb.



3.

Disinfection



Permitted: Alcohols, aldehydes, peroxide compounds (<5%) and quaternary compounds according to the disinfectants list of the Verbund für Angewandte Hygiene e.V. (Association for Applied Hygiene)

Recommended: Commercially available disinfection wipes (for surfaces). The surfaces of the silicone band and the buckle must be completely wetted and the movable part of the buckle should be moved over the entire length of the band. Attend contact time (usually 1 minute).

→ To prevent germs from being carried from patient to patient!