

Victoria RA631



Victoria radiation protection vest – Optimized fitted comfort and protection

The Victoria radiation protection vest has been specially developed for women to offer optimum protection and maximum comfort. As part of the Victoria collection, the vest combines modern functionality with an elegant design to meet the needs of women in interventional radiology.

Here are the key features:

- ✓ **Perfect Fit:** The vest fits ideally and is particularly comfortable to wear without the need for a full front overlap.
- ✓ **Compatible:** It can be easily combined with our radiation protection skirt (see costume RA631).
- ✓ **Comfortable and Practical:** Equipped with shoulder pads as a standard, the vest offers maximum comfort, especially when worn for long periods of time.
- ✓ **Certified Protection:** Type testing and certification in accordance with (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

Hint: Please note that sleeves cannot be attached to the Victoria vest due to the adapted armholes. Please combine the Victoria vest with the bolero RA637.

Protection/Lead Equivalent Standard Collection

Novalite® Plus or Leadfree 2.0



Front:

0.25 mm Pb, 0.35 mm Pb or 0.50 mm Pb

Back generally:
0.25 mm Pb

Protection/Lead Equivalent Levio Edition

Novalite® Plus or Leadfree 2.0



Total value Levio-Panel:

0.42 mm Pb or 0.50 mm Pb (outside 0.25 mm Pb)

Back generally:
0.25 mm Pb



Sizes*, recommended for:

| | | |
|-----|----------------------------------|----------------------------------|
| XS | Chest girth ≤ 84 cm | Waist girth ≤ 86 cm |
| S | Chest girth approx. 85 - 95 cm | Waist girth approx. 87 - 97 cm |
| M | Chest girth approx. 96 - 106 cm | Waist girth approx. 98 - 107 cm |
| L | Chest girth approx. 107 - 117 cm | Waist girth approx. 108 - 118 cm |
| XL | Chest girth approx. 118 - 128 cm | Waist girth approx. 119 - 129 cm |
| XXL | Chest girth approx. 129 - 139 cm | Waist girth approx. 130 - 140 cm |
| 3XL | Chest girth approx. 140 - 150 cm | Waist girth approx. 141 - 151 cm |

* For the vest, the larger value of the chest and waist girth determines the size.



ALL-ROUND PROTECTION: VEST ESPECIALLY FOR THE FEMALE BODY

VICTORIA RA631

Experience the difference with the Victoria radiation protection vest - the ideal combination of functionality, comfort and modern design, especially for Her needs in radiology.

Levio® 
Optimised Protection Concept

Innovative Design

I. Dedicated Fit for Women

The Victoria collection offers innovative X-ray protection that is specially adapted to the female anatomy. Adapted armholes protect the sensitive breast tissue without restricting freedom of movement.

II. Princess Seams These ensure a body-hugging, perfect fit and contribute to a comfortable and secure fit of the vest.

III. Stretch Insert in the Back The stretch insert at the back provides additional comfort and flexibility.

Levio® Edition

The vest RA631 is also available in our Levio Edition. Optimized weight-to-protection ratio: Instead of a uniform lead equivalent across the entire front area, Levio protects - from above the chest to below the gonads - with a higher lead equivalent.

Shoulder Relief System

Combined the vest with our shoulder relief system BRYGGA.



RFID chip included

Higher Protection – Levio Edition

The area with the reinforced protective material is visually separated/marked by a sewn-in band.

Optimised Armholes

An elastic band is incorporated into the armhole from size XL upwards to ensure that it fits perfectly in all sizes and thus provides the best possible protection for the breast tissue.

Dividing Seams in the Front

For a good fit to achieve full radiation protection

Stretch Insert

For a comfortable fit

Stretch Belt, Optionally

Additionally, a 10 cm wide, button-on stretch belt can be added to the vest.



Outer Material

Colours
(see front page)



Standard:
ComforTex® HPMF Hybrid

Sizes

from XS to 3XL
(see front page)



Optional

On request:
Stretch belt (external)
Breast pocket
Embroidery
Made-to-measure protection



Accessories/Replacement

Bolero
Thyroid/sternum protection



Stretch belt vest
(Adjustment via
snap fastener)