



Quanta System

# Youlaser Prime<sup>®</sup>

---

mixed technology

AESTHETICS & DERMATOLOGY

Quanta System

# Youlaser Prime<sup>®</sup>

---

mixed technology

## BRING YOUR SKIN TO ITS PRIME

Youlaser Prime ensures the already affirmed **high-quality standards**, offering **exclusive and refined treatment solutions** thanks to the reinvented Quanta Mixed Technology<sup>®</sup>.

## SYNERGY = PROGRESS

Youlaser Prime is the laser device you choose when you decide to upgrade your practice by **combining ablative and non-ablative emissions simultaneously**. This capability significantly **expands treatment possibilities**.

## ENHANCED FLEXIBILITY

This innovative laser device sets a new **standard of flexibility**, thanks to the **second generation of Quanta Mixed Technology<sup>®</sup>** and **three new exclusive treatments**: StarGaze, PrimeLift and IconicPeel.

## YOUR SKIN IN ITS PRIME

**Modularity and versatility** enable Youlaser Prime to achieve **balanced treatment results**, minimize downtime, and **prioritize consumers well-being**. This way, everyone's natural beauty can shine brightly.

## REACH THE EXCELLENCE

Thanks to its **significant heritage** and the **innovative technology** that characterizes it, the new Youlaser Prime is the ideal laser device **for all those operators who are seeking to extend their results**.

# Overview

## Two Combinable Wavelengths

CO<sub>2</sub> 10600 nm

Fiber Laser 1550 nm

## The Hyperpulse Power

Up to 150 W with CO<sub>2</sub> 10600 nm,  
for remarkable precision

## Quanta Mixed Technology®


The 2nd generation of Quanta Mixed Technology® allows for an enhanced combination of wavelengths through the selection of three mixing modes

## Three Exclusive Treatments


Exclusively available with Youlaser Prime, the IconicPeel, Stargaze, and PrimeLift treatments are designed to address each consumer's unique skin needs

## Multiple Handpieces Availability

New surgical handpieces, OptiScan Prime and Shelase Prime, have been designed to address a large variety of applications



FLEXIBILITY HAS NO BOUNDARIES



Synergy =  
Progress



# Innovation renews beauty

At Quanta System, we are dedicated to raising awareness about the importance of **embracing diversity and celebrating the unique beauty of every individual**.

**Youlaser Prime** fully reflects the demands of modern aesthetic medicine, representing **a high-level solution in beauty treatments** with **synergy as its driving force**.

This **innovative laser device** is the result of a profound commitment to **technological advancement** and is already establishing itself as **a concrete solution for addressing a wide range of skin needs and concerns**.



## MAIN TECHNOLOGICAL ADVANCEMENTS

ABLATIVE  
+ NON ABLATIVE  
COMBINATION

150W  
HYPERPULSE  
CO<sub>2</sub> 10600 nm

UP TO 18W  
WITH  
FL 1550 nm

2nd GEN  
QUANTA MIXED  
TECHNOLOGY

EXPERIENCE THE EXTRAORDINARY

# Refine your — practice

Continuously driven by **technological innovation**, Quanta System has developed an update to its proprietary **Quanta Mixed Technology®**, aiming to **enhance client comfort**.

This technology allows the **combination of ablative and non ablative laser wavelengths** in a single pulse, offering **a wide range of applications**.

Youlaser Prime's refined **Quanta Mixed Technology®** allows for the **mixed emission of two wavelengths from different sources**: **CO<sub>2</sub> 10600 nm** and **Fiber Laser 1550 nm**.


Enhancing precision over individual emissions, this advanced technology enables to **achieve a remarkable level of control**, offering three combination modes: **MIX mode**, **Hybrid mode**, and the original **MIX-Hybrid mode**.



YOULASER PRIME  
Quanta Mixed Technology®

CO<sub>2</sub> 10600 nm

FIBER LASER 1550 nm



THE RIGHT COMBINATION IS IN YOUR HANDS



Enhanced  
Flexibility

# Care for your — skin needs

Youlaser Prime<sup>®</sup>  
mixed technology

Quanta System's focus on client wellbeing has led to the design of **three exclusive treatments**, all of which can be performed only with Youlaser Prime and offer **a refined solution to many of today's most common skin concerns**.


These special procedures stand as a confirmation of the **technological innovation** behind this device, which builds on **Youlaser Family's decade-long legacy**, evolves in its Prime version, and unlocks **new treatment possibilities**.

**StarGaze**, **PrimeLift**, and **IconicPeel** provide three distinct, tailored solutions that **address the unique needs of consumers, helping them feel confident in their own skin and embrace the power of glowing up**.




#### STARGAZE

Complete periocular resurfacing



#### PRIMELIFT

Facial lifting effect induced by the Quanta Mixed Technology<sup>®</sup>, through the Hybrid mode




#### ICONICPEEL

Advanced peeling treatment

EXCLUSIVE SOLUTIONS FOR UNIQUE NEEDS





# Your skin in its — Prime

**Youlaser Prime** is dedicated to enhancing beauty and wellness by offering **sophisticated solutions that carefully meet customers' needs**, with the aim of celebrating and preserving **natural inner beauty**.

Quanta System's **profound commitment to continuous improvement** has led to a further upgrade of the proprietary **Quanta Mixed Technology®**, always placing **consumer's comfort at the heart of every treatment**.

Together, we can reimagine the meaning of **"beauty"**, reminding ourselves of the importance of **listening to our needs and practicing self-care**.

Elevate your skin to its prime.

**Skin Resurfacing/Rejuvenation**

---

**Surgical Treatments**

---

**Feminine Treatments**

---



# Treatments



## Skin Resurfacing/Rejuvenation

---

Combining two wavelengths simultaneously allows for a variety of **skin resurfacing and rejuvenation** procedures on **specific areas of face and neck**.

Thanks to controlled thermal diffusion, **scars, fine lines and wrinkles** can be treated efficiently, addressing the specific needs of each consumer.

**Rapid healing** is ensured due to the **minimal skin damage** that normally occurs during the procedure.



## Surgical Treatments

---

The Youlaser Family enables the delivery of a variety of surgical treatments, allowing users to perform procedures such as **blepharoplasty, rhinophyma, and actinic keratosis**.

**Multiple spot sizes** are available to optimize laser emission for the **removal of various lesions**, significantly **improving patients' skin health**.

To deliver these treatments, practitioners need to use **specific handpieces to treat skin tissues at deep levels**.





## Feminine Treatments

**Genitourinary syndrome of menopause (GSM)**, previously known as **atrophic vaginitis** or **vulvovaginal atrophy**, affects more than half of women in their **postmenopausal phase**. This condition, as well as **urinary incontinence** and **general vaginal laxity**, can have serious implications for **women's quality of life and sexual health**.

**Shelase** is an **innovative vaginal rejuvenation laser technology** designed by Quanta System to **reduce vaginal atrophy and laxity**, along with **enhance vulvar aesthetics**. Vaginal atrophy, characterized by the **loss of natural collagen thickness and structural integrity** in the vaginal tissue, can occur as a natural consequence of aging or childbirth.

Shelase utilizes a **dual laser technology** for an effective treatment that involves both **non-ablative and ablative resurfacing of the vaginal mucosa**. This approach results in **significant improvements in patients' quality of life** with minimal downtime, allowing for a rapid return to everyday activities.

This feminine treatment is performed to **stimulate contraction of mucosa and production of new collagen**, to regain the normal **vaginal elasticity**. Shelase procedures allow women to regain their confidence and take care of their **intimate health**.

# Handpieces

## OptiScan PRIME

Designed to reduce downtime, it provides both speed and adaptability, enhancing workflows

### Treatments Options

Skin Resurfacing  
Skin Rejuvenation  
StarGaze  
PrimeLift  
IconicPeel

INTEGRATED  
TOUCHSCREEN DISPLAY

QUICK  
CONNECTOR

CONTACT  
STAINLESS STEEL  
END-TIP

FINGER SWITCH



## Surgical Handpieces

Featuring innovative design, lightweight construction, and compact size

### Available Spot Sizes

150  $\mu\text{m}$   
200  $\mu\text{m}$   
400  $\mu\text{m}$   
1000  $\mu\text{m}$   
2000  $\mu\text{m}$



# Handpieces

## Shelase PRIME

Dedicated to women's intimate health, it delivers effective vaginal rejuvenation treatments

### Available Attachments

Speculum 90 mm

Speculum 120 mm

Vulvar Tip (for vulva rejuvenation and reshape)

Vaginal Ring (for using the vulvar tip without any speculum)









Reach the  
Excellence

# Challenge your — boundaries

Youlaser Prime has all the key features needed to **expand aesthetic treatment possibilities**, enabling operators to perform **high-level procedures** and improve their performances.

Thanks to its **remarkable potential** and **advanced emission control**, Youlaser Prime is **the perfect addition to any clinic seeking a significant turning-point** in their practice.

This **innovative laser device** is the result of a deep commitment to **technological advancement** and demonstrates to be a **concrete solution for addressing a wide range of skin needs and concerns**.




*CONTROLLED  
THERMAL  
DIFFUSION*

*EASY-TO-USE  
HANDPIECE  
DESIGN*

*SPECIFIC  
TREATMENT  
DATABASE*

A MILESTONE TO ACHIEVE SUCCESS

# Before/After



**STARGAZE**  
Complete periocular resurfacing  
*Courtesy of Quanta System*



**PRIMELIFT**  
Facial lifting effect induced  
by the Quanta Mixed  
Technology®, through the  
Hybrid mode  
*Courtesy of Quanta System*



**ICONICPEEL**  
Advanced peeling treatment  
*Courtesy of Quanta System*



# References

## **Youlaser Prime**

Al-Dhalimi, M. A., Al-Janabi, M. H. & Hussein, R. A. A. A. *The Use of a 1,540nm Fractional Erbium-Glass Laser in Treatment of Alopecia Areata*. *Lasers Surg Med* 51, 859–865 (2019).

Asri, E., Widayati, R. & Malik, D. *Effectiveness of combined 10,600 nm fractional CO<sub>2</sub> and 1540 nm erbium GaAs laser therapy on acne scar score alteration in patients with atrophic acne scars*. *Turk. J. Dermatol.* 13, 126 (2019).

Mezzana, P. *“Two wavelengths endovaginal laser system”: Clinical evaluation of a new device for mild SUI and vaginal atrophy treatment*. *Dermatol Ther* 33, e14445 (2020).

Mezzana, P., Valeriani, M. & Valeriani, R. *Combined fractional resurfacing (10600 nm/1540 nm): Tridimensional imaging evaluation of a new device for skin rejuvenation*. *J. Cosmet. Laser Ther.* 18, 397–402 (2016).

Paradisi, A., Ricci, F. & Sbano, P. *“Drug-resistant granuloma faciale”: treatment with carbon dioxide-GaAs laser*. *Dermatol Ther* 29, 317–319 (2016).

Paradisi, A., Ricci, F. & Sbano, P. *“Severe rhinophyma treated with double wavelengths carbon dioxide-GaAs laser: a case series”*. *Int J Dermatol* 60, e7–e8 (2021).

## **Shelase Prime**

Salinas Pena, J., Tameish, S., Guilarte Calzada, C., & Cavallé Busquets, P. (2025). *Efficacy of a Mixed Wavelength Laser for Vaginal Health in Postmenopausal Women: A Randomized Controlled Trial*. *International Journal of Women’s Health*, 17, 571–584. <https://doi.org/10.2147/IJWH.S486323>

# Advantages

## OUTSTANDING FLEXIBILITY

With reinvented Quanta Mixed Technology® and Hybrid/Mix emission mode, Youlaser Prime stands out for its flexibility

## ADVANCED CONTROL

The controlled combination of ablative and non-ablative emissions, delivers remarkable precision and effectiveness

## PRIME TREATMENT EXPERIENCE

Youlaser Prime responds to various skin needs of consumers, providing a refined treatment experience and concrete results

## REMARKABLE PERFORMANCES

With its upgraded features, Youlaser Prime represents a high-end solution in the aesthetic market

MORE INFO HERE



get in touch



[quantasystem.com](https://www.quantasystem.com)



This brochure is not intended for the U.S. market.  
Certain Intended Uses/Configurations/Models/Accessories are not cleared for U.S.

#### **DANGER - VISIBLE AND INVISIBLE LASER RADIATION**

Avoid eye or skin exposure to direct or scattered radiation | Laser product: Class 4 | Aiming beam: Class 3R

