

Hand-held magnifier luminaire

OPTICLUX Hand 10-2 UV D16005000-00807073

Derungs

MEDICAL LIGHTING



fitted with	LED - Energy efficiency category A+
color rendering index (CRI)	> 95
colour temperatur	6500 K
work equipment	electronic converter
mains lead	1.5 m Euro
lens; dimensions	Ø 160 mm
lens; diopters	3,5
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	0.6 kg
Power consumption	7 W
usage	wireless switch
class of protection	II
norms	EN 60601-1
special features	Battery, Ergonomic handle, scratch-proof plastic lens, Battery charger, Protective cover, continuously dimmable, switchable Wood-Light
order number	D16005000 - pure white

light distribution

